

L34 ANSWER 1 OF 4 BIOTECHDS COPYRIGHT 2005 THE THOMSON CORP. ON STN
DUPLICATE 1

AN 2002-11189 BIOTECHDS

TI **Detection of Helicobacter pylori
catalase activity in biological samples for diagnosis
of Helicobacter pylori infection;
antibody and enzyme immobilization for bacterium disease**

diagnosis

AU SUZUKI N; WAKASUGI M; MOCHIDA R; HIRATA H
PA WAKAMOTO PHARM CO LTD

PI WO 2002014541 21 Feb 2002

AI WO 2000-JP6894 11 Aug 2000

PRAI JP 2001-123221 20 Apr 2001

DT Patent

LA Japanese

OS WPI: 2002-227351 [28]

AB DERMENT ABSTRACT:

NOVELTY - **Diagnosis of Helicobacter pylori**
infection by measurement of catalase activity in digestive
tract or feces specifically due to **Helicobacter**
pylori catalase.

DETAILED DESCRIPTION - **Diagnosis of Helicobacter**
pylori infection, comprising measurement of catalase
activity in digestive tract or feces specifically due to
Helicobacter pylori catalase. An INDEPENDENT
CLAIM is included for kits for the diagnostic method.

BIOTECHNOLOGY - Preferred Method: The **Helicobacter**
pylori catalase in the sample is separated from other
catalases by ion-exchange or an immunochemical method.

USE - The method is useful for simple and effective
detection of infection by **Helicobacter pylori**

EXAMPLE - Monoclonal antibody specific to

Helicobacter pylori catalase was immobilized
on a multiwell plate. A sample of feces (3g) was suspended in
phosphate buffer (12 ml) and centrifuged at 4 degrees C/3000 rpm for 15
minutes. The supernatant (0.2 ml/well) was added to the plate. After one
hour at room temperature the plate was washed and 5 mM hydrogen peroxide
(0.2 ml/well) was added. After one hour at room temperature, oxygen
evolution in the wells was measured using a Chemetrics K-7512 apparatus.
Oxygen evolution of 8-12 parts per million (ppm) or more was shown by
feces samples positive for the **Helicobacter catalase**,
while samples negative for the **catalase** had oxygen evolution
below 1 ppm. (18 pages)